



CooperVision®
Live Brightly.®

Product Specifications Guide



clariti® 1 day family

clariti® 1 day



Great value and high breathability¹ in a convenient daily disposable contact lens.

Wearing Schedule
Daily disposable

Sphere Power
+ 8.00 to -10.00
(to +/- 6.00 in
0.25 steps)
(from +/- 6.00 in
0.50 steps)

Carton Size Options
30 or 90 pack

Design
Aspheric

Material / H₂O Content
somofilcon A / 56%

Base Curve 8.6mm
Diameter 14.1mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint No

clariti® 1 day toric



Great value and high breathability¹ in a convenient daily disposable contact lens.

Wearing Schedule
Daily disposable

Sphere Power
Plano to -6.00D
(0.25DS steps)

-6.50D to -9.00D
(0.50DS steps)

Carton Size Options
30 pack

Cylinder Power
-0.75
-1.25
-1.75

Material / H₂O Content
somofilcon A / 56%

Axis
10°, 20°, 60°,
70°, 80°, 90°,
100°, 110°,
120°, 160°,
170°, 180°

Base Curve 8.6mm
Diameter 14.3mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint No

Design
Aspheric

+0.25D to +4.00D
(0.25DS steps)

-0.75
-1.25
-1.75

20°, 70°,
90°, 110°,
160°, 180°

Plano to -6.00D
(0.25DS steps)

-2.25

10°, 20°,
90°, 160°,
170°, 180°

-6.50D to -9.00D
(0.50DS steps)

clariti® 1 day multifocal



Great value and high breathability¹ in the convenience of a daily disposable contact lens.

Wearing Schedule
Daily disposable

Sphere Power
+5.00D to -6.00D
(in 0.25 Steps)

Carton Size Options
30 pack

Add Power
Low:
Add up to
+2.25D

High:
+2.50D to
+3.00D

Material / H₂O Content
somofilcon A / 56%

Design
Centre near,
progressive
intermediate, and
peripheral distance

Base Curve 8.6mm
Diameter 14.1mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint No

Biofinity[®] family

Biofinity[®]



Wearing Schedule

Daily or 6 nights/7 days Extended

Sphere Power

+8.00D to -12.00D
(0.50D steps after +/-6.00D)

Carton Size Options

3 or 6 pack

Design

Aspheric

Material / H₂O Content

comfilcon A / 48%

Base Curve 8.6mm

Diameter 14.0mm

FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Extraordinary visual performance.
Premium all-day comfort.

Biofinity[®] toric



Wearing Schedule

Daily or up to 7 days Extended

Sphere Power

+8.00D to -10.00D
(0.50D steps after +/-6.00D)

Carton Size Options

3 or 6 pack

Cylinder Power

-0.75
-1.25
-1.75
-2.25

Material / H₂O Content

comfilcon A / 48%

Axis

10° to 180°
(10° steps)

Design

Optimised Toric Lens Geometry™

Base Curve 8.7mm

Diameter 14.5mm


FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Extraordinary visual performance.
Premium all-day comfort.

Biofinity[®] multifocal



Wearing Schedule

Daily Wear or Extended Wear up to 6 nights/7 days

Sphere Power

+6.00D to -8.00D
(0.50D steps after -6.00D)

Carton Size Options

3 or 6 pack

Add Power

+1.00
+1.50
+2.00
+2.50

Material / H₂O Content

comfilcon A / 48%

Design

D lens
N lens

Base Curve 8.6mm

Diameter 14.0mm

FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Extraordinary visual performance at every distance. Premium all-day comfort.

Avaira Vitality™ family

Avaira Vitality™



High performance.
Outstanding value.

Wearing Schedule

Daily

Sphere Power

-12.00D to +8.00D
(0.50D steps after
+/- 6.00D)

No Plano

Carton Size Options

3 or 6 pack

Design

Aspheric front surface

Material / H₂O Content

fanfilcon A / 55%

Base Curve 8.4mm
Diameter 14.2mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes

Avaira Vitality™ toric



High performance.
Outstanding value.

Wearing Schedule

Daily

Sphere Power

+0.25 to +6.00DS
(0.25D steps)

+6.50 to +8.00DS
(0.50D steps)

-0.25 to -6.00DS
(0.25D steps)

-6.50 to -10.00DS
(0.50D steps)

Carton Size Options

3 or 6 pack

Cylinder Power

-0.75
-1.25
-1.75
-2.25

Material / H₂O Content

fanfilcon A / 55%

Axis

10° to 180°
(10° steps)

Design

Optimised Toric
Lens Geometry™

Base Curve 8.5mm
Diameter 14.5mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes

Proclear® family

Proclear® 1 day



Clear vision², convenience and comfort, for contact lens wearers who experience eye dryness.

Wearing Schedule
Daily disposable

Carton Size Options
30 or 90 day pack

Material / H₂O Content
omafilcon A / 60%

Base Curve 8.7mm
Diameter 14.2mm
FDA Group 2
UV Blocking* No
Visibility Tint Yes

Sphere Power
+8.00D to -12.00D
(0.50D steps after -6.00D and +5.00D)

Design
Aspheric

Proclear® 1 day multifocal



Convenience and comfort for multifocal contact lens wearers who experience eye dryness.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
omafilcon A / 60%

Base Curve 8.7mm
Diameter 14.2mm
FDA Group 2
UV Blocking* No
Visibility Tint Yes

Sphere Power
+6.00D to -10.00D
(0.50D steps after -6.00D)

Add Power
Single power profile, up to +2.50

Design
Centre near aspheric

Proclear®



Clear vision³ and comfort for contact lens wearers who experience eye dryness.

Wearing Schedule
Daily

Carton Size Options
3 or 6 pack

Material / H₂O Content
omafilcon B / 62%

Base Curve 8.6mm
Diameter 14.2mm
FDA Group 2
UV Blocking* No
Visibility Tint Yes

Sphere Power
+20.00D to -20.00D (0.50D steps after +/- 6.00D)

Design
Sphere

Plano to -10.00D
(0.50D steps after +/- 6.00D)

Proclear® toric



Helps provide all-day comfort for contact lens wearers with astigmatism who experience eye dryness.

Wearing Schedule	Carton Size Options	Material / H ₂ O Content	Base Curve
Daily	3 or 6 pack	omafilcon B / 62%	8.4, 8.8mm
Sphere Power	Cylinder Power	Axis	Diameter 14.4mm
+6.25D to +10.00D (0.50D steps after +6.50D)	-0.75 -1.25 -1.75 -2.25	5° to 180° (5° steps)	FDA Group 2
-8.50D to -10.00D (0.50D steps)			UV Blocking* No
			Visibility Tint Yes
+10.00D to -10.00D (0.50D steps after +/- 6.50D)	-2.25 to -5.75 (0.50D steps)	5° to 180° (5° steps)	
+6.00 to -6.25D (0.25D steps) -6.50 to -8.00D (0.50D steps)	-2.75 to -5.75DC (0.50D steps)	5° to 180° (5° steps)	
+6.25, +6.50 to +10.00D (0.50D steps)	-0.75 to -5.75DC (0.50D steps)	5° to 180° (5° steps)	
-8.50 to -10.00D (0.50D steps)	-0.75 to -5.75DC (0.50D steps)	5° to 180° (5° steps)	

Proclear® multifocal



Helps provide all-day comfort for multifocal contact lens wearers who experience eye dryness.

Wearing Schedule
Daily

Carton Size Options
3 or 6 pack

Material / H₂O Content
omafilcon B / 62%

Base Curve 8.7mm

Diameter 14.4mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power

+6.00D to -8.00D
(0.50D steps after
-6.50D)

Add Power

+1.00
+1.50
+2.00
+2.50D

Design

D lens
N lens

+6.25D to +20.00D
-8.50D to -20.00D
(0.50D steps after
+/-6.50D)

+1.00
+1.50
+2.00
+2.50

+20.00D to
-20.00D
(0.50D steps after
+/-6.50D)

+3.00
+3.50
+4.00

+20.00D to
-20.00D
(0.50D steps after
+/-6.50D)

+1.00
+1.50
+2.00
+2.50
+3.00
+3.50
+4.00

Proclear® multifocal toric



Helps provide all-day comfort for multifocal contact lens wearers with astigmatism who experience eye dryness.

Wearing Schedule
Daily

Carton Size Options
3 or 6 pack

Material / H₂O Content
omafilcon B / 62%

Base Curve 8.4, 8.8mm

Diameter 14.4mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power

+20.00D to
-20.00D (0.50D
steps after +/-
6.50D)

Cylinder Power

-0.75 to -5.75
(0.50 steps)

Axis

5° to 180°
(5° steps)

Design

D lens
N lens

Biomedics® family

Biomedics® 55 Evolution



Wearing Schedule	Carton Size Options	Material / H₂O Content	Base Curve	8.8mm
Daily	6 pack	ocufilcon D / 55%	Diameter	14.2mm
Sphere Power	Design		FDA Group	4
-0.50D to -10.00D (0.50D steps after -6.00D)	Aspheric		UV Blocking*	Yes
			Visibility Tint	Yes

Designed for comfort and performance at an affordable price.


Biomedics® toric



Wearing Schedule	Carton Size Options	Material / H₂O Content	Base Curve	8.7mm
Daily	3 or 6 pack	ocufilcon D / 55%	Diameter	14.5mm
Sphere Power	Cylinder Power	Axis	FDA Group	4
+6.00D to -9.00D (0.50 steps after -6.00D and +5.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	UV Blocking*	Yes
			Visibility Tint	Yes

The affordable contact lens for wearers with astigmatism.

Biomedics® 1 day extra



Wearing Schedule	Carton Size Options	Material / H₂O Content	Base Curve	8.6mm
Daily	30 pack	ocufilcon D / 55%	Diameter	14.2mm
Sphere Power	Design	Axis	FDA Group	4
-0.50D to -10.00D (0.50D steps after -6.00D)	Sphere	10° to 180° (10° steps)	UV Blocking*	Yes
			Visibility Tint	Yes

Greater comfort, preferred by wearers.⁴

* UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.

¹ Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygenation During Contact Lens Wear. Optom Vis Sci 2005; 82(6):467-472.

² CVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

³ Due to aspheric optics.

⁴ Vs Biomedics 1 day. Clinical Performance of the Biomedics 1 day Extra vs Biomedics 1 day, 2011, data on file.